


Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

INSTRUCTIONS FOR USE

This manual **MUST** be read **BEFORE** using this product



Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

Contents

1.	Introduction.....	3
2.	Symbols	3
3.	Product Labelling.....	4
4.	Warnings and Precautions for Use	7
5.	Product Overview	9
6.	Product Features	10
7.	Product Description.....	11
8.	Installation	12
9.	Product Operation.....	13
	9.1 Control Panel	13
	9.2 Control Unit Operation	13
	9.3 Mattress operation.....	15
10.	Cleaning and Disinfection.....	19
	10.1 Cleaning and Disinfection Protocol: Control Unit.....	19
	10.2 Cleaning and Disinfection Protocol: Mattress	20
11.	Storage	21
12.	Troubleshooting	22
13.	Care and Preventative Maintenance.....	23
14.	Warranty	24
15.	Disposal	24
16.	Technical Specifications	25
	16.1 Control Unit.....	25
	16.2 Mattress and Top Cover	25
	16.3 Transport and Operating Conditions.....	26
	16.4 Safety Standards.....	26
	16.5 Electromagnetic Compatibility	26
17.	Contact Information	31



Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

1. Introduction

This document outlines important information and instructions for use (IFU) pertaining to the safe and effective use of the product. Read all instructions carefully before using the product. Store the IFU in a designated area, where it is always easily accessible. If unsure, consult a medical professional regarding the correct use of the product. For further product related information, contact Winncare PAC Ltd directly; see the “Contact Information” section of this document.

2. Symbols










Symbols and advisory notices are used in this document to help safe and optimal operation of the product. See information below for definitions of the symbols.

 WARNING	<p>Warning: Safety warning. Failure to obey and understand could lead to injury to yourself or others, and in some circumstances death.</p>
 CAUTION	<p>Caution to highlight potential hazards that, if failed to follow, could lead to damage or failure in parts or all of the system and equipment.</p>
<div style="border: 1px solid black; padding: 5px; display: inline-block;">NOTE</div>	<p>Note: Important information users should be aware of for correct use of the equipment.</p>














Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

3. Product Labelling









The labels shown are for illustrative purposes only – some symbols on your product may differ from the examples shown.

	<p>Instructions for Use Read the Instructions for use before use</p>
	<p>Type BF Applied Part Applied Part: The parts of the device that come into physical contact with the user/occupant in order for it to carry out its intended function. Type BF: Applied parts which are electrically isolated from earth and other parts of the medical equipment - Complying with specific requirements for protection against electric shock to IEC 60601-1.</p>
	<p>W.E.E.E Label Waste Electrical and Electronic Equipment.</p>
	<p>Class II electrical device The user/occupant is protected by at least two layers of insulation between the current carrying parts (e.g. mains cable) – If damage is noticed to the control unit or mains cable assembly turn off at the mains supply and contact your provider or Winnocare PAC Ltd. immediately</p>
IP21	Protected from touch by fingers and objects greater than 12mm. Protected from condensation.
	CE marking indicating conformity with European Community harmonized legislation. Figures indicate Notified Body supervision.
	UK marking indicating conformity with UK Medical Device Regulations 2002 (SI 2002 No 8, as amended)
	Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745.
	Serial number
	Reference number

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)


	Batch code
	Name and address of the manufacturer
	Date of manufacture
	Temperature limitation to indicate the temperature limitation for the product during usage
	Unique device identifier
	Authorized Representative in the European Community
	Disinfect by wiping the surface using a hypochlorite solution diluted 1000 ppm.
	Machine wash up to 95°C.
	Tumble dry on a low setting
	Do not use harsh abrasives or Phenol cleaners
	Do not iron
	Ensure system is dry before storing, use and reuse.
	Do not place heavy objects on surface of cover other than the patient

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

	<p>Do not use when damp, ensure surface is dry before use</p>
	<p>Do not fold. Roll pack the system</p>
	<p>Do not use sharp objects</p>
	<p>Max Patient weight defines the maximum total load of the patient kg (lb)</p>
	<p>Safe Working load (SWL) is the maximum combined weight of the patient and any equipment that the mattress can safely support.</p>
	<p>Foot end</p>
	<p>Resistant to ignition</p>
	<p>Recycling</p>

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)


4. Warnings and Precautions for Use


WARNING


Do not use device control unit in oxygen rich environment or near flammable gases. **RISK OF FIRE AND BURN INJURY.**


WARNING

Do not use device with a damaged power cable. **RISK OF ELECTROCUTION AND FIRE.**


WARNING

Ensure appropriate cable management. Avoid operating the device with loose or severely taught cables. **RISK OF TRIP AND FALL INJURY.**


WARNING


Do not open or repair the control unit whilst it is in use or connected to mains power supply. **RISK OF ELECTRIC SHOCK.**


WARNING


Do not use the device as a repositioning tool. **RISK OF PRESSURE INJURY.**


WARNING


CPR connector must be accessible at all times. **RISK OF SERIOUS INJURY OR DEATH.**


WARNING


Ensure the device is assembled and operated as intended. **RISK OF PRESSURE INJURY.**


WARNING


Do not cover the control unit with blankets and other items. **RISK OF FIRE.**


WARNING

Do not spray liquid on the control unit whilst it is connected to mains power. **RISK OF ELECTRICAL BURNS.**










WARNING

Do not expose any parts of the device to a naked flame. Do not smoke. **RISK OF FIRE AND PROPERTY DAMAGE.**


WARNING

Ensure the patient is manually repositioned at frequent intervals. **RISK OF PRESSURE INJURY.**

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

 CAUTION Complete a risk assessment if any accessory is used with the device.	 CAUTION Ensure there are no additional layers between the surface of the mattress and the patient. Device performance may be affected. Complete a risk assessment if in doubt.
 CAUTION To ensure optimal function of the device, use suitably trained personnel for servicing and repair. Use original parts only.	 CAUTION Ensure the device is suitable for the patient. Complete a risk assessment if in doubt. Consult a medical professional if in doubt.
 CAUTION Ensure the device is plugged into mains power supply for optimal function.	 CAUTION Ensure the mains cable does not become compressed, trapped or damaged by the bed frame or other equipment.
 CAUTION Do not use device alongside hot water bottles or electric blankets. Device performance may be affected.	 CAUTION Complete a risk assessment when using device with incontinence products.

NOTE	Use a CE marked extension cable if device power cable cannot reach wall socket. If in use, do not overload.
-------------	---


Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

5. Product Overview

This product is made of the following components:


Control Unit

1. Control Panel
2. On/Off Button
3. Female Air Connector Port
4. Air Filter
5. Fuse Holder
6. Pad
7. Hooks



Mattress

1. Top Cover
2. Air Cells
3. Air tube assembly leading to Male Air Connector (not shown)
4. CPR Connector at end of air tube set (not shown)
5. Sloped heel section



Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

6. Product Features

HER-COMBO-777-CLP: Moray Combo Hybrid Mattress System CLP	
Size (mm) L/W/H	1980 x 880 x 150
Maximum Patient Weight	< 230 kg (36 stone)
Safe working load	230 kg (36 stone)
Other Features	<ul style="list-style-type: none"> • One in two cell cycle design (AB cell pattern) • 9 cells: 9 fully alternating cells • 4" Foam in cell design • Nylon PU air cell construction • Sloped visco memory foam profiled heel zone • Welded, multi-stretch, waterproof, vapour permeable cover with antibacterial and antifungal properties • Machine washable cover up to 95°C • CPR connector (on mattress hose) for rapid deflation • 11.5 mins cycle time • Pressure range: 12-28 mmHg • Audible low-pressure alert • Semi-auto pressure adjustment • Auto Firm function with auto-return • Alert mute • Control panel lock • Constant Low Pressure mode (CLP)

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

7. Product Description

Intended Use	To provide pressure relief and aid in the prevention and management of pressure related injuries as part of a standard package of care.
Target Population	Typical adults with limited mobility, undergoing some medical supervision and monitoring. Individuals assessed as “at risk” and up to “very high risk” of pressure damage and/or with existing tissue damage, as determined by a combination of clinical judgment and validated assessment tools.
Contraindications	Patients below the minimum or maximum user weight listed for the associated device. Cervical or skeletal traction. Unstable skeletal fractures. Unstable spinal injury.
Users	Caregivers, laypersons and/or medical professionals.
Warranty	2 years subject to regular maintenance and servicing.
Reusable	Devices are re-usable and must be cleaned in between each patient use.
Maintenance or calibration	Perform regular mattress audits to check for fluid ingress and strike-through on mattress top cover. The system should be serviced once a year, as a minimum.
Accessories	Devices are not sold with accessories.
Risk Assessment	It is the responsibility of the end user/care provider to carry out the necessary risk assessment to ensure the patient’s safety. This should be carried out before using the mattress system. A risk assessment should include, but is not limited to: <ul style="list-style-type: none"> • Product combinations (bed frame, mattress, side rails etc.) • Extent of tissue damage (if any) • Entrapment • Patient falls • Small adults (and children) • Patients with learning difficulties • Patients with atypical anatomy • Unauthorised people with access to the controls • Use with other medical accessories e.g. incontinence products

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

8. Installation

To install the device, follow the instructions below:

1. Carefully open the packaging.
2. Check the product for any signs of damage. Do not use if damaged and contact your provider or WinnCare PAC Ltd.
3. Place the mattress on top of the bed frame with the top cover facing upwards and the air tube set at the foot end of the bed.
4. Attach the mattress to the bed frame by securing with the adjustable securing straps (if applicable).

NOTE


On profiling beds, it is essential that adjustable straps are secured around the movable sections of the bed frame, otherwise the mattress may be damaged.

5. Ensure the CPR connector (on mattress hose) is easily accessible at all times.
6. Using the hooks on the back of the control unit, hang the unit over the frame/board at the foot end of the bed. If there is no foot frame/board lay the unit on the floor, under the bed with the front control panel facing upwards.
7. Attach the air tube set using the male air connector on the mattress to the female air connector port on the control unit, ensuring the air tubing is not kinked or trapped between parts of the bed frame/other equipment.
8. Plug the mains cable into a suitable mains supply and switch on the control unit.
9. The mattress will start to inflate when the control unit is switched on.
10. Once fully inflated, adjust the straps that attach the mattress to the bed frame, ensuring the mattress is held securely in place (if applicable)
11. Cover the mattress loosely with a sheet, ensuring it does not interfere with cell alternation.







Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

9. Product Operation



9.1 Control Panel




9.2 Control Unit Operation

 Power	Power On / Off	 SWL	Safe Working Load: 230kg
 Auto Firm	<p>All cells will inflate to the firmest setting of 28 mmHg (with no dynamic alternation – also known as max inflate) & the digital display will show “- - -”. The control unit will revert back to the previous operating mode after 30 minutes.</p>	 Lock / Unlock	<ul style="list-style-type: none"> • Push the button for 2-3 seconds to lock / unlock the panel. • The LED will flash continuously: <ul style="list-style-type: none"> ○ when the mattress has not yet reached pressure during alternating cycles. ○ If the mattress has an air leak, or the pressure has dropped lower than 6 mmHg. • When the panel is unlocked, it will automatically lock after a few minutes. • The lock/unlock function will not operate if the system has a pressure failure.
 Mute	This button will silence the audible alert when pressed		Press to increase or decrease the pressure mode. The current mode will be indicated on the LCD display.

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

 Engineer	<p>The Engineer LED will flash upon initialisation (the digital display will also flash). Once the system has reached pressure the LED will extinguish (digital display will now be a solid light). If Initialisation fails (approx. 50 minutes), the LED will continue to flash and the audible alert will sound. The Engineer LED will also flash if the pressure has dropped below 6mmHg for 6 minutes. (Please see the 'Incomplete inflation/low pressure' section of the Troubleshooting guide)</p>	Digital Display & Pressure Settings	
			The digital display will flash when initialising (along with the engineer light). Once the system has reached pressure and is ready to use the engineer light will extinguish and the digital display light will be solid (no flashing).
		1	CLP: Constant Low Pressure at 16 mmHg with no dynamic alternation. There is no auto-return, the control unit will remain in CLP mode until changed.
		2	Soft: Alternating mode at pressure of 12 mmHg
		3	Medium: Alternating mode at pressure of 18 mmHg
		4	Firm: Alternating mode at pressure of 24 mmHg
5	Total Support: Alternating mode at pressure of 28 mmHg. Can be used as seating mode to increase the pressure when backrest is profiled.		



CAUTION

Ensure there is no accidental deactivation of the system. Function lock does not lock the on/off switch.

NOTE	Auto Firm will automatically reset to alternating mode after 30 minutes.
-------------	---

NOTE	The mute setting will automatically cancel after 15 minutes. The audible signal will re-sound.
-------------	---

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

9.3 Mattress operation


1. Turn on the power to the control unit. The control unit will start to inflate the mattress.
2. The lock/unlock LED will continually flash whilst the system reaches pressure. Once the system reaches pressure (after approx.. 30-40 minutes) the Ready For Use LED will illuminate.
3. The system will default to the last setting used. This must be manually changed to suit the individual patient by using the comfort control buttons.
4. If the setting selected is not providing a comfortable and supportive pressure, then the pressure can be increased or decreased accordingly to ensure a more comfortable and supportive pressure using the Comfort Control. Switch off the function lock and adjust the pressure to provide a comfortable pressure level for the patient.
5. Using clinical judgement and with continuous monitoring of the patient for up to **72** hours, increase or decrease the pressure levels to suit the patients comfort levels.

Auto Firm Mode


For patient safety, the mattress should always be used in alternating mode. Auto firm mode may be selected for short periods to provide a firm base for clinical/nursing needs.

When auto firm mode is selected, all air cells inflate to the pressure set, creating a static surface.

If, after 30 minutes, the control unit is still set to auto firm mode, it will automatically return to alternating mode. This is for patient safety, to ensure they are not left on a constantly inflated surface.




Alternating Mode



Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)


Constant Low Pressure (CLP) Mode

When the system is in CLP mode, the mattress remains at a constant pressure of 16mmHg with no dynamic alternation. This helps to envelope & immerse the user by increasing the surface area over which the user is supported. There is no auto-return and the control unit will remain in CLP mode until changed. CLP is sometimes referred to as static mode.




Total Support (TS) Mode

When Total Support is selected, the control unit increases the pressure to 28mmHg whilst alternating. This setting can also be used as a seating mode to increase the pressure in the system when the backrest is profiled.



Profiled



NOTE

The mattress can be used on a profiling bed where the backrest is profiled to an angle of 65°. Pressure settings may need to be increased. Use clinical judgement to optimise pressure relief.



Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

Semi-Auto Pressure Adjustment

Once the pressure level has been set, the control unit monitors the mattress pressure and maintains it at the set level. If the pressure falls below this level, the control unit will automatically speed up inflation of the mattress until the correct pressure is met. If the control unit is unable to maintain the set pressure, an audible alert will sound, and the low-pressure indicator light will flash. If this occurs, refer to the 'Troubleshooting' section.

CPR Function



In an emergency rapid deflation of the mattress may be required. The CPR connector is located at the foot end of the mattress on the end of the air tube set. Disconnect from control unit for CPR.



To re-inflate reconnect the CPR connector back into the control unit. The mattress will start to inflate. Wait for optimal pressure to be reached before using the mattress.

Transporting the Mattress

If the mattress is disconnected from the power supply so it can be moved, or in the event of a mains power failure, carry out the following procedure to maintain mattress inflation:



1. Disconnect the air tube set from the control unit by squeezing the two tabs (A) on the male connector and pull away from the control unit (B).
2. Seal using the cap marked "Transport" which is attached via a cord to the male connector.
3. Switch off the control unit.
4. Disconnect from the power supply.
5. The mattress can now be moved.

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

NOTE

In the event of a power fail, the mattress will remain inflated for up to 24 hours. The mattress should be returned to the mains supply as soon as possible. If not plugged into the mains device performance will be affected.

NOTE

Do not drag the mattress, always carry it.

NOTE

Complete instructions above quickly to minimize air loss



WARNING

Do not remove the mattress from the bed frame if the occupant is still on the mattress. RISK OF FALL.









WARNING

When not connected to mains power supply, alternating mode will not available. RISK OF LIMITED PRESSURE RELIEF.

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

10. Cleaning and Disinfection

Cleaning and disinfection of the device is mandatory. As a minimum, cleaning and disinfection should be performed before the device is taken out of storage, between each patient, at regular intervals whilst in use and before being placed into storage. Consult local practice guidelines for more details on cleaning and disinfection of reusable medical devices. Follow the instructions below to achieve a minimum level of cleaning and disinfection for safe use of the device.

 <p>WARNING</p> <p>Wear appropriate personal protective equipment (PPE) when cleaning the mattress or control unit. RISK OF SKIN IRRITATION</p>	 <p>CAUTION</p> <p>Do not immerse or soak the control unit. Do not spray any cleaning solution onto the control unit.</p>
 <p>CAUTION</p> <p>Do not attempt to clean the device whilst its connected to mains power.</p>	 <p>CAUTION</p> <p>Do not use Phenol-based solutions or abrasive compounds.</p>
 <p>CAUTION</p> <p>Do not use top cover if strike-through or damage is suspected.</p>	 <p>CAUTION</p> <p>Do not autoclave.</p>

10.1 Cleaning and Disinfection Protocol: Control Unit

1. Visually check the product for external damage – do not use if damage is found.
2. Place the control unit on a work surface and using a clean, soft, non-abrasive cloth, wipe the outside of the case with a prepared sodium hypochlorite solution (recommended 1,000 ppm).
3. The control unit should be cleaned by starting with the cleanest areas and systematically moving to the dirtiest areas. Extra care should be taken around areas where excess dirt or dust may gather.
4. Change the cloth if it becomes dirty.
5. Once clean, wipe down with a fresh, clean, soft, non-abrasive cloth, moistened with clean water to remove detergent residue.
6. Dry off with a paper towel. Always allow the surfaces to dry thoroughly before putting back into use.

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

10.2 Cleaning and Disinfection Protocol: Mattress

NOTE

Before attempting to clean the mattress, the top cover should be checked for physical signs of damage that may lead to strike-through (ingress of fluid through cover). Staining to the underside of the top cover is a sign of strike-through.

1. The mattress should be regularly checked for damage or tears. Replace if damaged.
2. Wipe down with a clean cloth moistened with a mild detergent and diluted in warm water (40°C).
3. Rinse with cold clean water using a clean, soft, non-abrasive cloth and allow to fully dry before use.

Disinfection

1. Unzip the top cover from the mattress.
2. The top cover can be machine washed up to 95°C and tumble dried on a cool setting.
3. Carefully clean all air cells with (1,000 ppm) prepared solution of sodium hypochlorite and allow to dry completely. In extreme circumstances 10,000 ppm can be used, wipe with cold water to finish (Frequent cleaning with a high concentration disinfectant solution (i.e. 10,000 ppm available chlorine) may reduce the life span of the system).
4. Re-assemble the mattress.
5. Ensure the mattress is completely dry before either storing or reusing.


NOTE

Frequent or prolonged exposure to higher concentration disinfectant solutions may prematurely age the fabric cover of the mattress.

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)


11. Storage

Follow the instructions below to prepare the product for storage.




CAUTION

Ensure device is cleaned and disinfected prior to storage.




CAUTION

Do not fold, crease or stack mattress.



CAUTION

Do not stack control units with other medical equipment.



CAUTION

Do not stack the control units when in storage.



CAUTION


Do not store whilst inflated.

1. Detach the control unit from the mattress.
2. Store in a sealed polythene bag to protect from dirt, debris, fluids etc. with a suitable identification tag.
3. Store the control unit in a sealed polythene bag to protect from dirt, debris, fluids etc. with a suitable identification tag.

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

12. Troubleshooting

Table below provides a guide to the product’s alerts.



WARNING

**Do not open or repair the control unit whilst it is in use or connected to mains power supply.
RISK OF ELECTRIC SHOCK.**

NOTE	If mains cable or plug is visibly damaged turn off power supply at the mains and contact your approved services engineer.
-------------	---

Problem	Actions
Power Failure	<ol style="list-style-type: none"> 1. Turn off the control unit to silence the alarm and unplug from the mains supply <div style="border: 1px solid black; padding: 5px; margin: 10px 0;"> <p style="text-align: center;">NOTE</p> <p>If the ‘power failure’ indicator activates the mute button will not silence the audible signal. To silence, the system must be turned off by pressing the power button.</p> </div> <ol style="list-style-type: none"> 2. Check the mains socket is working - plug in a device that is known to work. 3. Plug the control unit back into the wall socket. 4. Turn on the control unit. If the control unit still fails to operate: 5. Turn off the control unit at the wall & replace plug fuse. 6. Turn on the control unit. If the control unit still fails to operate: 7. Replace control unit fuses. For fuse types see ‘Technical Specifications’ section. 8. Turn on the control unit. 9. If control unit still fails to operate, turn off at the mains and contact your approved service provider.
Incomplete inflation/low pressure	<ol style="list-style-type: none"> 1. Ensure the mattress air connector is properly connected to the control unit, is not constricted in any way and has no kinks. 2. Turn the unit off and then on again to clear the indicator. If the ‘low pressure’ indicator continues to illuminate. 3. Remove the top cover and ensure there is no air leakage within the mattress – cells, tubing and connectors. 4. Turn the unit off and then on again to clear the indicator.

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

Problem	Actions
	5. If a low-pressure indicator is still evident turn off at the mains and contact your approved service provider
Alternating mode failure	<ol style="list-style-type: none"> 1. Turn off the control unit. 2. Disconnect the male air connector to reduce cell pressure. 3. Reconnect air connector. 4. Turn on the control unit. 5. If alternating mode is still inoperable turn off at the mains and contact your approved service provider
Patient is bottoming out	<ol style="list-style-type: none"> 1. Ensure the patient is suited to the rating of the mattress. 2. Ensure the patient is centrally positioned on the mattress. 3. Increase the pressure setting – refer to ‘Mattress Operation’ section. 4. If the patient is still bottoming out refer to ‘incomplete inflation’ above

13. Care and Preventative Maintenance

The expected service life of this product is 6 years, subject to appropriate servicing and use in accordance with these instructions. Winnicare PAC Ltd. recommends annual servicing of this product as a minimum. For optimal performance of the device, more frequent visual and operational inspections are encouraged wherever possible. Contact Winnicare PAC Ltd to arrange your annual service. Failure to do so may invalidate product warranty.

NOTE	Always disconnect the control unit from the mains power supply prior to performing any maintenance procedures (when viable).
NOTE	No modification of this equipment is allowed. Use original parts only.
NOTE	The mattress system should be vacated by the patient before any maintenance or inspection takes place. If this is not possible due to the patient’s mobility, care should be taken for the service engineer not to make contact with the patient when working on electrical items.

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

14. Warranty

This product is covered by manufacturer's warranty as part of the General Terms and Conditions of Business. Some warranty periods may differ – refer to the product features section of this document for the exact warranty period pertaining to your product.

Any warranty claims during the warranty period must be investigated by Winncare PAC Ltd, where return of the original product maybe required. A warranty claim is successful if the product is faulty due to a manufacturing defect. This warranty does not cover any other damage including but not limited to; misuse, natural wear and tear, lack of maintenance, accidental damage and unauthorized modifications.

15. Disposal

Should the product reach the end of its use and no longer be repaired, ensure that it is disposed of in accordance with local W.E.E.E. (Waste Electrical and Electronic Equipment) policies. Alternatively, contact Winncare PAC Ltd to arrange for collection. The metal and plastic components used in both the mattress and control unit should be separated and recycled – consult local recycling practices for further information.



Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

16. Technical Specifications

16.1 Control Unit

Control Unit: PUMP1031	
Dimensions (mm) H x W x D	140 x 260 x 95
Weight (kg)	1.35
Cycle time (min)	11.5
Air Output (L/min)	>5
Power cord length (m)	4.5
Noise Level	≤45dB
Supply Rating	220 ~ 240V, 50Hz, 8W
Fuse Rating	Mains Plug – 3A Control Unit – F1A, 250V
Mains Plug	Type G/BS1363/A
Electrical classification	Electrical shock protection: Type BF Applied Part: Liquid ingress protection: IP21 Not AP or APG equipment*
<i>*Not suitable for use in the presence of flammable aesthetic mixtures with air, oxygen or nitrous oxide.</i>	

16.1.1 Internal Cell Pressure

According to the weight data manually selected to suit the individual patient, the table below highlights the corresponding internal cell pressure.

CLP (1)	16mmHg	Soft (2)	12mmHg	Medium (3)	18mmHg
Firm (4)	24mmHg	Total Support (5)	28mmHg	Auto Firm	28mmHg

16.2 Mattress and Top Cover

Mattress Base Unit: HER-COMBO-777-NP and Top Cover: COVER6022 & COVER6083	
Number of cells	9
Cell Material	Nylon PU
Cell depth (inches)	4" Foam in Cell
Base Cover Material	Nylon PVC
Internal Foam	Base Foam: 38/39 kg U-Core Heel Foam: 50 kg Visco Cell Foam: 38 kg & 40/41 kg & 50 kg Visco
Weight (kg)	15
Emergency	CPR connector
Top Cover Material	PU and Polyester Mix

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

16.3 Transport and Operating Conditions

Transport and storage conditions	Temperature: -25°C to +70°C Humidity: < 93% max, non-condensing
Operational conditions	Temperature: +5°C to +40°C Humidity: 15% - 93%, non-condensing Atmospheric Pressure: 700hPa to 1060hPa Operating Altitude: ≤ 2000m Pollution: Degree 2 UV: Intended for indoor use only

16.4 Safety Standards

BS EN 60601-1:2006+A13:2024
IEC 60601-1-11
IEC 60601-1-2
BS EN 61000
IEC 61000-3-3
IEC 61000-3-2
IEC 61000-4-2
IEC 61000-4-4
IEC 61000-4-5
IEC 61000-4-11
IEC 61000-4-8
IEC 61000-4-6
IEC 61000-4-3

16.5 Electromagnetic Compatibility

The control unit has been designed to meet the EMC requirements of BS EN 61000. This standard defines the levels of immunity to electromagnetic interferences as well as maximum levels of electromagnetic emissions for medical devices. They are in place to provide reasonable protection against dangerous interference in a medical or residential environment.

Immunity to electromagnetic interference - this refers to the levels of electromagnetic interference that the control unit can withstand from nearby sources radiating radio frequency (RF) energy (e.g. from mobile phones, network devices etc).

Electromagnetic emissions - this refers to the levels of RF energy the control unit emits.

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

The immunity levels are set out in the following manufacturers guidance. If these levels are exceeded, then the system may not operate correctly or stop operating. It is important therefore to try to ascertain the source of the interference by turning nearby equipment off. There are simple measures that can be taken to correct the problem:

- Remove or relocate the interfering equipment,
- Increase the separation distance between the control unit and the interfering equipment.

The RF emissions are set out in the following manufacturers guidance. The control unit generates very low RF energy, however interference to sensitive equipment is still possible. If interference to radio/tv reception and/or other equipment is suspected, turning the control unit off and on can determine if this is the case. There are simple measures that can be taken to correct the problem:

- Relocate the receiving antenna,
- Increase the separation distance between the control unit and affected equipment.

Due to the increasing number of wireless devices, such as laptops and mobile phones, it is important that the system is installed following the manufacturer’s guidance to ensure continued and reliable operation.


Requirements according to BS EN 61000 and BS EN 60601-1:2006+A13:2024: Guidance and manufacturer’s declaration – electromagnetic emissions		
Emission test	Compliance	Electromagnetic environment – guidance
RF emission CISPR 11	Group 1	The control unit uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	Dynamic mattress system is suitable for use in all establishments, including domestic establishments and those directly connected to the public, low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Requirements according to BS EN 61000 and BS EN 60601-1:2006+A13:2024: Guidance and manufacturer’s declaration – electromagnetic immunity			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

Requirements according to BS EN 61000 and BS EN 60601-1:2006+A13:2024: Guidance and manufacturer's declaration – electromagnetic immunity			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s)	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT† (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
<i>†UT is the a.c. mains voltage prior to application of the test level.</i>			
Conducted RF IEC 61000- 4-6 Radiated RF IEC 61000- 4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 Vrms 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the control unit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance. d = 1.2√P d = 1.2√P 80 MHz to 800 MHz d = 2.3√P 800 MHz to 2.5 GHz

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

Requirements according to BS EN 61000 and BS EN 60601-1:2006+A13:2024: Guidance and manufacturer's declaration – electromagnetic immunity			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
			<p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey*, should be less than the compliance level in each frequency range**.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>
<p>* Field strengths from fixed transmitters cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Dynamic mattress system is used exceeds the applicable RF compliance level above, Dynamic mattress system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the system.</p> <p>** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

NOTE	<p>At 80 MHz and 800 MHz, the higher frequency range applies. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>
-------------	---

The dynamic mattress system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer/user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the dynamic mattress system as recommended below, according to the maximum output power of the communications equipment.

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

Recommended separation distances between portable and mobile RF communications equipment and the control unit			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m) Electromagnetic environment – guidance		
	150 KHZ TO 80 MHZ D = 1.2VP	80 MHZ TO 800 MHZ D = 1.2VP	800 MHZ TO 2.5 GHZ D = 2.3VP
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE

At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Moray Combo Hybrid Mattress System CLP (HER-COMBO-777-CLP)

17. Contact Information



WINNCARE PAC Ltd.
Unit 2 Interchange Way,
Oakenshaw
Bradford,
BD12 7AZ
United Kingdom

Customer Service
Tel +44 0800 193 6030
Email: askus@winncare.uk
Servicing Email: askservice@winncare.uk
www.winncare.uk



WINNCARE France
4, Le Pas du Chateau,
85670 Saint-Paul-Mont-Penit
France
Tel. +33 (0)4 66 02 15 15

www.winncare.fr



NOTE

Inform the competent authority if you believe or have reason to believe that the device presents a serious risk or that it has been tampered with.
All serious incidents which are related to the device must be notified to the manufacturer and the competent authority of the member state in which the user and/or patient resides.