

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

### INSTRUCTIONS FOR USE

This manual **MUST** be read **BEFORE** using this product



# Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

## Contents

1.	Introduction.....	3
2.	Symbols .....	3
3.	Product Labelling.....	4
4.	Warnings and Precautions for Use .....	7
5.	Product Overview .....	9
6.	Product Features .....	10
7.	Product Description.....	11
8.	Installation .....	12
9.	Product Operation.....	13
	9.1 Control Panel .....	13
	9.2 Control Unit Operation .....	13
	9.3 Mattress operation .....	15
10.	Cleaning and Disinfection.....	19
	10.1 Cleaning and Disinfection Protocol: Control Unit.....	19
	10.2 Cleaning and Disinfection Protocol: Mattress .....	20
11.	Storage.....	21
12.	Troubleshooting .....	22
13.	Care and Preventative Maintenance.....	23
14.	Warranty .....	23
15.	Disposal .....	24
16.	Technical Specifications .....	25
	16.1 Control Unit.....	25
	16.2 Mattress and Top Cover .....	25
	16.3 Transport and Operating Conditions .....	26
	16.4 Safety Standards.....	26
	16.5 Electromagnetic Compatibility .....	26
17.	Contact Information .....	31



# Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

## 1. Introduction

This document outlines important information and instructions for use (IFU) pertaining to the safe and effective use of the product. Read all instructions carefully before using the product. Store the IFU in a designated area, where it is always easily accessible. If unsure, consult a medical professional regarding the correct use of the product. For further product-related information contact Winncare PAC Ltd directly; see the “Contact Information” section of this document.

## 2. Symbols










Symbols and advisory notices are used in this document to help safe and optimal operation of the product. See information below for definitions of the symbols.

 <b>WARNING</b>	<p><b>Warning:</b> Safety warning. Failure to obey and understand could lead to injury to yourself or others, and in some circumstances death.</p>
 <b>CAUTION</b>	<p><b>Caution</b> to highlight potential hazards that, if failed to follow, could lead to damage or failure in parts or all of the system and equipment.</p>
<div style="border: 1px solid black; padding: 5px; display: inline-block;">NOTE</div>	<p><b>Note:</b> Important information users should be aware of for correct use of the equipment.</p>

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

### 3. Product Labelling









The labels shown are for illustrative purposes only – some symbols on your product may differ from the examples shown.

	<p><b>Instructions for Use</b> Read the Instructions for use before use</p>
	<p><b>Type BF Applied Part</b> Applied Part: The parts of the device that come into physical contact with the user/occupant in order for it to carry out its intended function. Type BF: Applied parts which are electrically isolated from earth and other parts of the medical equipment - Complying with specific requirements for protection against electric shock to IEC 60601-1.</p>
	<p><b>W.E.E.E Label</b> Waste Electrical and Electronic Equipment.</p>
	<p><b>Class II electrical device</b> The user/occupant is protected by at least two layers of insulation between the current carrying parts (e.g. mains cable) – If damage is noticed to the control unit or mains cable assembly turn off at the mains supply and contact your provider or Winnicare UK Ltd. Immediately.</p>
<p>IP21</p>	<p>Protected from touch by fingers and objects greater than 12mm. Protected from condensation.</p>
	<p>CE marking indicating conformity with European Community harmonized legislation. Figures indicate Notified Body supervision.</p>
	<p>UK marking indicating conformity with UK Medical Device Regulations 2002 (SI 2002 No 8, as amended).</p>
	<p>Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745.</p>
	<p>Serial number</p>
	<p>Reference number</p>

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)












	Batch code
	Name and address of the manufacturer
	Date of manufacture
	Temperature limitation to indicate the temperature limitation for the product during usage
	Unique device identifier
	Authorized Representative in the European Community
	Disinfect by wiping the surface using a hypochlorite solution diluted 1000 ppm.
	Machine wash up to 95°C.
	Tumble dry on a low setting
	Do not use harsh abrasives or Phenol cleaners
	Do not iron
	Ensure system is dry before storing, use and reuse.
	Do not place heavy objects on surface of cover other than the patient

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

	<p>Do not use when damp, ensure surface is dry before use</p>
	<p>Do not fold. Roll pack the system</p>
	<p>Do not use sharp objects</p>
	<p>Max Patient weight defines the maximum total load of the patient kg (lb)</p>
	<p>Safe Working load (SWL) is the maximum combined weight of the patient and any equipment that the mattress can safely support.</p>
	<p>Foot end</p>
	<p>Resistant to ignition</p>
	<p>Recycling</p>

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

### 4. Warnings and Precautions for Use

 <p><b>WARNING</b></p> <p>Do not use device control unit in oxygen rich environment or near flammable gases. <b>RISK OF FIRE AND BURN INJURY.</b></p>	 <p><b>WARNING</b></p> <p>CPR dial must be accessible at all times. <b>RISK OF SERIOUS INJURY.</b></p>
 <p><b>WARNING</b></p> <p>Do not use device with a damaged power cable. <b>RISK OF ELECTROCUTION AND FIRE.</b></p>	 <p><b>WARNING</b></p> <p>Ensure the device is assembled and operated as intended. <b>RISK OF PRESSURE INJURY.</b></p>
 <p><b>WARNING</b></p> <p>Ensure appropriate cable management. Avoid operating the device with loose or severely taught cables. <b>RISK OF TRIP AND FALL INJURY.</b></p>	 <p><b>WARNING</b></p> <p>Do not cover the control unit with blankets and other items. <b>RISK OF FIRE.</b></p>
 <p><b>WARNING</b></p> <p>Do not open or repair the control unit whilst it is in use or connected to mains power supply. <b>RISK OF ELECTRIC SHOCK.</b></p>	 <p><b>WARNING</b></p> <p>Do not spray liquid on the control unit whilst it is connected to mains power. <b>RISK OF ELECTRICAL BURNS.</b></p>
 <p><b>WARNING</b></p> <p>Do not use the device as a repositioning tool. <b>RISK OF PRESSURE INJURY.</b></p>	 <p><b>WARNING</b></p> <p>Do not expose any parts of the device to a naked flame. Do not smoke. <b>RISK OF FIRE AND PROPERTY DAMAGE.</b></p>
 <p><b>WARNING</b></p> <p>Ensure the patient is manually repositioned at frequent intervals. <b>RISK OF PRESSURE INJURY.</b></p>	

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)



**CAUTION**

Complete a risk assessment if any accessory is used with the device.



**CAUTION**

Ensure there are no additional layers between the surface of the mattress and the patient. Device performance may be affected. Complete a risk assessment if in doubt.



**CAUTION**

Ensure the device is suitable for the patient. Complete a risk assessment if in doubt. Consult a medical professional if in doubt.



**CAUTION**

To ensure optimal function of the device, use suitably trained personnel for servicing and repair. Use original parts only.



**CAUTION**

Ensure the mains cable does not become compressed, trapped or damaged by the bed frame or other equipment.



**CAUTION**

Ensure the device is plugged into mains power supply for optimal function.



**CAUTION**

Ensure the mains cable does not become compressed, trapped or damaged by the bed frame or other equipment.



**CAUTION**

Do not use device alongside hot water bottles or electric blankets. Device performance may be affected.



**CAUTION**

Complete a risk assessment when using device with incontinence products.

**NOTE**

Use a CE marked extension cable if device power cable cannot reach wall socket. If in use, do not overload.

# Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

## 5. Product Overview

This product is made of the following components:

Pure Air 8 AR Mattress System (PUR-8-AR)	
Control Unit	CU-PUR-4
Mattress Base Unit	PUR-8-AR-M
Top Cover	SMP515

### Control Unit

1. Control Panel
2. On/Off switch
3. Mains Power Cable
4. Female Air Connector Port
5. Air Filter
6. Fuse Holders
7. Pad
8. Hooks



### Mattress

1. Top Cover
2. Air Cells
3. Air Tube Set with Male Air Connector
4. Easy-turn CPR Dial
5. Base Cover
6. Securing Straps



## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

### 6. Product Features

Pure Air 8 AR Mattress System (PUR-8-AR)	
Size (mm) L/W/H	2000 x 900 x 203
Maximum Patient Weight	250 kg (40 stone)
Safe working load	250 kg (40 stone)
Other Features	<ul style="list-style-type: none"> <li>• One in two cell-cycle design giving optimum therapy</li> <li>• 20 fully alternating air cells</li> <li>• 8" depth, split cell design including static base</li> <li>• PU air cell construction</li> <li>• Welded, PU multi-stretch waterproof and vapour permeable cover</li> <li>• Machine washable cover up to 95°C</li> <li>• Easy turn CPR dial for rapid deflation</li> <li>• Optional 10/15 mins cycle time</li> <li>• Pressure range: 20-45 mmHg</li> <li>• Silent running pump at optimum support pressures</li> <li>• Audible low-pressure alert</li> <li>• Audible high-pressure alert</li> <li>• Automatic pressure adjustment</li> <li>• Optional manual pressure adjustment</li> <li>• Static function with auto-return</li> <li>• CLP function</li> <li>• Cycle fault audible alert</li> <li>• Cycle time control</li> <li>• Power fail alert</li> <li>• Alert mute</li> <li>• Control panel lock</li> </ul>

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

### 7. Product Description

<b>Intended Use</b>	To provide pressure relief and aid in the prevention and management of pressure related injuries as part of a standard package of care.
<b>Target Population</b>	Typical adults with limited mobility, undergoing some medical supervision and monitoring. Individuals assessed as “at risk” and up to “very high risk” of pressure damage and/or with existing tissue damage as determined by a combination of clinical judgment and validated assessment tools.
<b>Contraindications</b>	Patients below the minimum or maximum user weight listed for the associated device. Cervical or skeletal traction. Unstable skeletal fractures. Unstable spinal injury.
<b>Users</b>	Caregivers, laypersons and / or medical professionals.
<b>Warranty</b>	3 years subject to regular maintenance and servicing.
<b>Reusable</b>	Devices are re-usable but must be cleaned in between each patient use.
<b>Maintenance or calibration</b>	Perform regular mattress audits to check for fluid ingress and strike-through on mattress top cover. The system should be serviced once a year, as a minimum.
<b>Accessories</b>	Devices are not sold with accessories.
<b>Risk Assessment</b>	It is the responsibility of the end user/care provider to carry out the necessary risk assessment to ensure the patient’s safety. This should be carried out before using the mattress system. A risk assessment should include, but is not limited to: <ul style="list-style-type: none"> <li>• Product combinations (bed frame, mattress, side rails etc.)</li> <li>• Extent of tissue damage (if any)</li> <li>• Entrapment</li> <li>• Patient falls</li> <li>• Small adults (and children)</li> <li>• Patients with learning difficulties</li> <li>• Patients with atypical anatomy</li> <li>• Unauthorized people with access to the controls</li> <li>• Use with other medical accessories e.g. incontinence products</li> </ul>

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

### 8. Installation

To install the device, follow the instructions below:

1. Carefully open the packaging.
2. Check the product for any signs of damage. Do not use if damaged and contact your provider or Winncare PAC Ltd.
3. Place the mattress on top of the bed frame with the top cover facing upwards and the male air connector at the foot end of the bed.
4. Attach the mattress to the bed frame by securing with the adjustable securing straps.

**NOTE**

On profiling beds, it is essential that adjustable straps are secured around the movable sections of the bed frame, otherwise the mattress may be damaged.

5. Rotate the CPR dial (located at the foot end of the mattress) to a vertical or horizontal, closed position (as below):




6. Using the hooks on the back of the control unit, hang the unit over the frame/board at the foot end of the bed. If there is no foot frame/board lay the unit on the floor, under the bed, with the front control panel facing upwards.
7. Attach the air tube set using the male air connector on the mattress to the female air connector port on the control unit/pump, ensuring the air tubing is not kinked or trapped between parts of the bed frame/other equipment.
8. Plug the mains cable into a suitable mains supply and switch on the control unit.
9. The mattress will start to inflate and will be completely inflated within 30-45 minutes.
10. Once fully inflated, adjust the straps that attach the mattress to the bed frame, ensuring the mattress is held in place securely.
11. Cover the mattress loosely with a sheet, ensuring it does not interfere with cell alternation.






# Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)


## 9. Product Operation

### 9.1 Control Panel



### 9.2 Control Unit Operation





	<p><b>1. Automatic Pressure Adjustment</b> An amber indicator illuminates when in automatic pressure adjustment mode (default mode). The amber indicator will switch off when in manual pressure adjustment mode.</p>
	<p><b>2. Power Failure</b> An amber indicator flashes and an audible alert will sound if a power failure occurs.</p>
 	<p><b>3. Low Pressure/Cycle Fault</b> An amber indicator illuminates when setting the pressure initially. An amber indicator flashes and an audible alert sounds if the pressure becomes unacceptably low during operation.</p> <p>An amber indicator illuminates and an audible alert sounds if alternation failure occurs.</p>
	<p><b>4. Function Lock</b> The control unit will automatically lock out all functionality 2 minutes after a function change. To unlock the control unit, the 'lock' button is pressed for 2 seconds. To re-engage the lock, the button can either be pressed again for 2 seconds or the user can wait for the automatic lock to re-engage. When the system is locked an amber indicator illuminates.</p>









**CAUTION**

**Ensure there is no accidental deactivation of the system. Function lock does not lock the on/off switch.**

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

  	<p><b>5. Cycle Time</b> Press the button to change between cycle times. Refer to the label on the back of the control unit and top cover label for your product code.</p> <p>15 minutes</p> <p>10 minutes</p> <p>The green indicator next to the icons will be illuminated to show which one is selected.</p>
	<p><b>6. Static Mode</b> Selecting 'Static' mode fully inflates all cells with no dynamic alternation (also known as max inflate) – amber light illuminated.</p>


<b>NOTE</b>	Static mode will automatically reset to alternating mode after 30 minutes.
-------------	--

	<p><b>7. Mode Select/Comfort Adjustment</b> Press and hold to switch between automatic &amp; manual pressure adjustment modes. When in manual pressure adjustment mode, press to cycle through the pressure settings.</p> <p>In alternating mode there are 4 available pressure settings (20 mmHg, 25 mmHg, 35 mmHg, 45 mmHg ±2 mmHg).</p> <p>In CLP mode there is 1 available pressure setting: 15mmHg. The green light illuminates to indicate which of the settings is operational.</p>
	<p><b>8. Operation Modes – Alternating &amp; CLP Modes</b> Refer to the CE label on the back of the control unit for product code.</p> <p><b>Alternating Mode:</b> Selecting 'Alternating' mode  inflates and deflates the cells in sequence over the cycle time selected – green light illuminated.</p> <p><b>CLP (Constant Low Pressure) Mode:</b> Selecting 'CLP' mode inflates all  cells to deliver continuous low pressure increasing the area of support – amber light illuminated.</p> <p>Alternating Sequence (PUR-8-AR)</p> <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">  </div> <div style="text-align: center;">  </div> </div>

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

--	--

**NOTE** In CLP mode pressure remains at a constant 15mmHg

	<p><b>9. Alert Mute</b> To mute an audible alert, press the button. The amber indicator will illuminate. Re-press the button to reset the alert.</p>
---	--

**NOTE** The mute setting will automatically cancel after 15 minutes. The audible signal will re-sound.

**NOTE** If the 'power failure' indicator activates the mute button will not silence the audible signal. To silence, the system must be turned off by pressing the power button.

### 9.3 Mattress operation

1. Turn on the power on the control unit. The pump starts to inflate the mattress.
2. The mattress will inflate in default 'Alternating' automatic pressure adjustment mode. Four indicators on the pressure dial will illuminate. The amber low pressure/cycle fault indicator will illuminate as inflation commences.
3. Green indicators will illuminate to show the cycle time. The control unit defaults to a 10-minute cycle time, 25 mmHg pressure range and alternating mode.
- 4.

<b>Default mode:</b>	
Operation Mode:	Alternating
Alternating Cycle Time:	10 minutes
Pressure Setting:	25 mmHg
<b>Complete settings:</b>	
Operation Mode	Alternating
Alternating Cycle Time:	10 minutes 15 minutes
Pressure Setting:	20 mmHg (Stage 1) 25 mmHg (Stage 2) 35 mmHg (Stage 3) 45 mmHg (Stage 4)

4. After 2 minutes the amber function lock indicator will be illuminated.
5. Once optimum pressure is reached (about 30-45 minutes) the control unit will beep twice,

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

and the amber low pressure/cycle fault indicator will switch off. Once a patient is positioned on the mattress, the control unit will enter automatic pressure adjustment mode, and the pressure setting will be automatically determined by the patient’s position and morphology.

6. If automatic pressure adjustment mode is not required, switch off the function lock, switch to manual pressure adjustment mode, and adjust the pressure to provide a suitable and comfortable pressure level for the patient.
7. Using clinical judgement and continuous monitoring of the patient for up to 72 hours, increase or decrease the pressure levels to suit the patients comfort levels.

**NOTE**

The mattress can be used on a profiling bed where the backrest is profiled to an angle of 65°. Pressure settings may need to be increased. Use clinical judgement to optimise pressure relief.

### Automatic Pressure Adjustment

Once the pressure level has been set, the control unit monitors the mattress pressure and maintains it at the set level. If the pressure falls below this level, the control unit will automatically increase the inflation of the mattress until the correct pressure is achieved. If the control unit is unable to maintain the set pressure, an audible alert will sound, and the low-pressure indicator light will flash. If this occurs refer to the ‘Troubleshooting’ section.

### CPR Function

In an emergency rapid deflation of the mattress may be required. The CPR dial is located at the foot end of the mattress. Rotate the CPR dial to the open position and the entire system will rapidly deflate.



To re-inflate turn the CPR dial to the closed position. The mattress will start to inflate. Wait for optimal pressure to be reached before using the mattress.

### CLP (Continuous Low Pressure) Mode

When CLP mode is selected, all air cells inflate to deliver continuous low pressure increasing the area of support, with the aim of reducing peak pressure at the interface via immersion and envelopment.

**NOTE**

For optimum performance in CLP mode, it is recommended that the surface be utilised in a recumbent position.

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

### Static Mode


For patient safety, the mattress should always be used in alternating mode. Static mode (max inflate) may be selected for short periods to provide a firm base for clinical/nursing needs.

When static mode is selected, all air cells inflate to the pressure that is set, creating a static surface. If, after 30 minutes, the control unit is still set to static mode it will automatically return to alternating mode. This is for patient safety, to ensure they are not left on a constantly inflated surface.

### Transporting the Mattress

If the mattress is disconnected from the power supply so it can be moved, or in the event of a mains power failure, carry out the following procedure to maintain mattress inflation:

1. Disconnect the air tube set with male connector from the power unit by squeezing the two tabs (A) and pulling away from the control unit (B).



2. Seal using the cap marked "Transport", which for safety is attached to the male connector.
3. Switch off the control unit.
4. Disconnect from the power supply.
5. The mattress can now be moved.

**NOTE**

In the event of a power fail, the mattress will remain inflated for up to 24 hours. The mattress should be returned to the mains supply as soon as possible. If not plugged into the mains device performance will be affected.

**NOTE**

Do not drag the mattress, always carry it.

**NOTE**

Complete instructions above quickly to minimize air loss

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)



**WARNING**

**Do not remove the mattress from the bed frame if the occupant is still on the mattress. RISK OF FALL.**










**WARNING**

**When not connected to mains power supply, alternating mode will not available. RISK OF LIMITED PRESSURE RELIEF.**

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

### 10. Cleaning and Disinfection

Cleaning and disinfection of the product is mandatory. As a minimum, cleaning and disinfection should be performed before the product is taken out of storage, between each patient, at regular intervals whilst in use and before being placed into storage. Consult local practice guidelines for more details on cleaning and disinfection of reusable medical devices. Follow the instructions below to achieve a minimum level of cleaning and disinfection for safe use of the device.

 <p><b>WARNING</b></p> <p><b>Wear appropriate personal protective equipment (PPE) when cleaning the mattress or control unit. RISK OF SKIN IRRITATION</b></p>	 <p><b>CAUTION</b></p> <p><b>Do not immerse or soak the control unit. Do not spray any cleaning solution onto the control unit.</b></p>
 <p><b>CAUTION</b></p> <p><b>Do not attempt to clean the device whilst its connected to mains power.</b></p>	 <p><b>CAUTION</b></p> <p><b>Do not use phenol-based cleaning solutions, solvents, neat bleach or abrasive products to clean the casing as this may cause damage.</b></p>
 <p><b>CAUTION</b></p> <p><b>Do not use top cover if strike-through or damage is suspected.</b></p>	 <p><b>CAUTION</b></p> <p><b>Do not use Phenol-based solutions or abrasive compounds.</b></p>
 <p><b>CAUTION</b></p> <p><b>Do not autoclave.</b></p>	

#### 10.1 Cleaning and Disinfection Protocol: Control Unit

1. Visually check the product for external damage – do not use if damage is found.
2. Place the pump on a work surface and, using a soft, clean, non-abrasive cloth, wipe the outside of the case with a prepared sodium hypochlorite solution (recommended 1,000 ppm).
3. The control unit should be cleaned by starting with the cleanest areas and systematically moving to the dirtiest areas. Extra care should be taken in areas where excess dirt or dust may gather.
4. Change the cloth if it becomes dirty.

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

5. Once clean, wipe down with a fresh, clean, soft, non-abrasive cloth moistened with clean water to remove detergent residue.
6. Dry off with a paper towel. Always allow the surfaces to dry thoroughly before reuse.

### 10.2 Cleaning and Disinfection Protocol: Mattress

**NOTE**

Before attempting to clean the mattress, the top cover should be checked for physical signs of damage that may lead to strike-through (ingress of fluid through cover). Staining to the underside of the top cover is a sign of strike-through.

1. The mattress should be regularly checked for damage or tears. Replace if damaged.
2. Wipe down with a clean cloth moistened with a mild detergent and diluted in warm water (40°C).
3. Rinse with cold, clean water and a clean cloth and allow to fully dry before use.

#### Disinfection

1. Unzip the top cover from the mattress.
2. The top cover can be machine washed up to 95°C and tumble dried on a cool setting.
3. Unsnap the air cells from the mattress base on both sides.
4. Carefully clean with (1,000 ppm) prepared solution of sodium hypochlorite and allow to dry completely. In extreme circumstances, 10,000 ppm can be used, wipe with cold water to finish (Frequent cleaning with a high concentration disinfectant solution (i.e. 10,000 ppm available chlorine) may reduce the life span of the system).
5. Make sure to disconnect all the air cells and spray the cleaning solution on all sides, including the connecting tubes and hoses.
6. Reassemble the mattress.
7. Ensure the mattress is completely dry before either storing or reusing.






**NOTE**

Frequent or prolonged exposure to higher concentration disinfectant solutions may prematurely age the fabric cover of the mattress.

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

### 11. Storage

Follow the instructions below to prepare the product for storage.


 <p><b>CAUTION</b></p> <p><b>Ensure device is cleaned and disinfected prior to storage.</b></p>	 <p><b>CAUTION</b></p> <p><b>Do not fold, crease or stack mattress.</b></p>
 <p><b>CAUTION</b></p> <p><b>Do not stack the control units when in storage.</b></p>	 <p><b>CAUTION</b></p> <p><b>Do not store whilst inflated.</b></p>
 <p><b>CAUTION</b></p> <p><b>Do not stack control units with other medical equipment.</b></p>	

1. Detach the control unit from the mattress.
2. Rotate the CPR dial until it is open.
3. Ensure there is no air trapped in the cells.
4. Lay the mattress out flat and roll the mattress from the foot end towards the head end.
5. Store in a sealed polythene bag to protect from dirt, debris, fluids, etc., with a suitable identification tag.
6. Store the control unit in a sealed polythene bag to protect from dirt, debris, fluids, etc., with a suitable identification tag.
7. If taking a device out of storage, unfold the mattress and allow it to lie unfolded for several minutes. Allow the product to acclimatise to the operating conditions.

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

### 12. Troubleshooting

Table below provides a guide to the product’s alarms.



**WARNING**

**Do not open or repair the control unit whilst it is in use or connected to mains power supply. RISK OF ELECTRIC SHOCK.**

<b>NOTE</b>	If mains cable or plug is visibly damaged turn off power supply at the mains and contact your approved services engineer.
-------------	---

<b>Problem</b>	<b>Actions</b>		
<b>Power Failure</b>	<ol style="list-style-type: none"> <li>1. Turn off the control unit to silence the alarm and unplug from the mains supply.</li> </ol> <div style="border: 1px solid black; padding: 5px; margin: 10px 0;"> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 15%; padding: 5px; text-align: center;"><b>NOTE</b></td> <td style="padding: 5px;">If the ‘power failure’ indicator activates the mute button will not silence the audible signal. To silence, the system must be turned off by pressing the power button.</td> </tr> </table> </div> <ol style="list-style-type: none"> <li>2. Check the mains socket is working - plug in a device that is known to work.</li> <li>3. Plug the control unit back into the wall socket.</li> <li>4. Turn on the control unit. If the control unit still fails to operate:</li> <li>5. Turn off the control unit at the wall &amp; replace the plug fuse.</li> <li>6. Turn on the control unit. If the control unit still fails to operate:</li> <li>7. Replace control unit fuses. For fuse types, see the ‘Technical Specifications’ section.</li> <li>8. Turn on the control unit. If the control unit still fails to operate, turn off at the mains and contact your approved service provider.</li> </ol>	<b>NOTE</b>	If the ‘power failure’ indicator activates the mute button will not silence the audible signal. To silence, the system must be turned off by pressing the power button.
<b>NOTE</b>	If the ‘power failure’ indicator activates the mute button will not silence the audible signal. To silence, the system must be turned off by pressing the power button.		
<b>Incomplete inflation/low pressure</b>	<ol style="list-style-type: none"> <li>1. Ensure the mattress air connector is properly connected to the control unit, is not constricted in any way and has no kinks.</li> <li>2. Ensure the CPR dial is closed and no air is leaking.</li> <li>3. Turn the unit off and then on again to clear the indicator. If the ‘low pressure’ indicator continues to illuminate:</li> <li>4. Remove the top cover and ensure there is no air leakage within the mattress cells, tubing and connectors.</li> <li>5. Turn the unit off and then on again to clear the indicator.</li> </ol>		

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

	If a low-pressure indicator is still evident, turn off at the mains and contact your approved service provider
<b>Alternating mode failure</b>	<ol style="list-style-type: none"> <li>1. Turn off the control unit.</li> <li>2. Disconnect the male air connector to reduce cell pressure.</li> <li>3. Reconnect the air connector.</li> <li>4. Turn on the control unit.</li> <li>5. If alternating mode is still inoperable, turn off at the mains and contact your approved service provider</li> </ol>
<b>Patient is bottoming out</b>	<ol style="list-style-type: none"> <li>1. Ensure the patient is suited to the rating of the mattress.</li> <li>2. Ensure the patient is centrally positioned on the mattress.</li> <li>3. Increase the pressure setting – refer to the ‘Mattress Operation’ section.</li> <li>4. If the patient is still bottoming out, refer to ‘incomplete inflation’ above</li> </ol>

### 13. Care and Preventative Maintenance

The expected service life of this product is 6 years, subject to appropriate servicing and use following these instructions. Winnicare PAC Ltd. recommends annual servicing of this product as a minimum. For optimal performance of the device, more frequent visual and operational inspections are encouraged wherever possible. Contact Winnicare PAC Ltd to arrange your annual service. Failure to do so may invalidate the product warranty.

<b>NOTE</b>	Always disconnect the control unit from the mains power supply prior to performing any maintenance procedures (when viable).
<b>NOTE</b>	No modification of this equipment is allowed. Use original parts only.
<b>NOTE</b>	The mattress system should be vacated by the patient before any maintenance or inspection takes place. If this is not possible due to the patient’s mobility, care should be taken for the service engineer not to make contact with the patient when working on electrical items.

### 14. Warranty



This product is covered by a manufacturer’s warranty as part of the General Terms and Conditions of Business. Some warranty periods may differ – refer to the product features section of this document for the exact warranty period of your product.

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

Any warranty claims during the warranty period must be investigated by Winncare PAC Ltd, where return of the original product may be required. A warranty claim is successful if the product is faulty due to a manufacturing defect. This warranty does not cover any other damage, including but not limited to: misuse, natural wear and tear, lack of maintenance, accidental damage and unauthorised modifications.

### 15. Disposal

Should the product reach the end of its use and may no longer be repaired, ensure that it is disposed of following local W.E.E.E. (Waste Electrical and Electronic Equipment) policies. Alternatively, contact Winncare PAC Ltd to arrange for collection. The metal and plastic components used in both the mattress and control unit should be separated, and recycled – consult local recycling practices for further information.



## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

### 16. Technical Specifications

#### 16.1 Control Unit

<b>Control Unit: CU-PUR-4</b>	
<b>Dimensions (mm) H x W x D</b>	205 x 280 x 115.8
<b>Weight (kg)</b>	2.6
<b>Cycle time (min)</b>	10 & 15 (AB)
<b>Air Output (L/min)</b>	6.5
<b>Power cord length (m)</b>	4.5
<b>Noise Level</b>	<40dB(A)
<b>Supply Rating</b>	230V, 50Hz, 12W
<b>Fuse Rating</b>	Mains Plug – 5A Control Unit - T1A, 250VAC
<b>Mains Plug</b>	Type G/BS1363/A
<b>Electrical classification</b>	Electrical shock protection: Class II Type BF Applied Part: Mattress Liquid ingress protection: IP21 Not AP or APG equipment*
<i>*Not suitable for use in the presence of flammable aesthetic mixtures with air, oxygen or nitrous oxide.</i>	

#### 16.2 Mattress and Top Cover

<b>Mattress Base Unit: PUR-8-AR-M and Top Cover: SMP515</b>	
<b>Number of cells</b>	20
<b>Cell Material</b>	TPU
<b>Cell depth (inches)</b>	8 (+ 4 inch static air base)
<b>Base Material</b>	Nylon PVC
<b>Weight (kg)</b>	10
<b>Emergency</b>	CPR dial
<b>Top Cover Material</b>	PU and Polyester mix

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

### 16.3 Transport and Operating Conditions

<b>Transport and storage conditions</b>	Temperature: -25°C to +70°C Humidity: < 93% max, non-condensing
<b>Operational conditions</b>	Temperature: +5°C to +40°C Humidity: 15% - 93%, non-condensing Atmospheric Pressure: 700hPa to 1060hPa Operating Altitude: ≤ 2000m Pollution: Degree 2 UV: Intended for indoor use only

### 16.4 Safety Standards

BS EN 60601-1:2006+A13:2024
IEC 60601-1-11
IEC 60601-1-2
BS EN 61000
IEC 61000-3-3
IEC 61000-3-2
IEC 61000-4-2
IEC 61000-4-4
IEC 61000-4-5
IEC 61000-4-11
IEC 61000-4-8
IEC 61000-4-6
IEC 61000-4-3

### 16.5 Electromagnetic Compatibility

The control unit has been designed to meet the EMC requirements of BS EN 61000. This standard defines the levels of immunity to electromagnetic interference as well as maximum levels of electromagnetic emissions for medical devices. They are in place to provide reasonable protection against dangerous interference in a medical or residential environment.

*Immunity to electromagnetic interference* - this refers to the levels of electromagnetic interference that the control unit can withstand from nearby sources radiating radio frequency (RF) energy (e.g. from mobile phones, network devices etc).

*Electromagnetic emissions* - this refers to the levels of RF energy the control unit emits.

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

The immunity levels are set out in the following manufacturer's guidance. If these levels are exceeded, then the system may not operate correctly or stop operating. It is important, therefore to try to ascertain the source of the interference by turning nearby equipment off. There are simple measures that can be taken to correct the problem:

- Remove or relocate the interfering equipment,
- Increase the separation distance between the control unit and the interfering equipment.

The RF emissions are set out in the following manufacturer's guidance. The control unit generates very low RF energy, however, interference to sensitive equipment is still possible. If interference to radio/tv reception and/or other equipment is suspected, turning the control unit off and on can determine if this is the case. There are simple measures that can be taken to correct the problem:

- Relocate the receiving antenna,
- Increase the separation distance between the control unit and affected equipment.

Due to the increasing number of wireless devices, such as laptops and mobile phones, it is important that the system is installed following the manufacturer's guidance to ensure continued and reliable operation.


Requirements according to BS EN 61000 and BS EN 60601-1:2006+A13:2024: Guidance and manufacturer's declaration – electromagnetic emissions		
Emission test	Compliance	Electromagnetic environment – guidance
RF emission CISPR 11	Group 1	The control unit uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	Dynamic Mattress System is suitable for use in all establishments, including domestic establishments and those directly connected to the public, low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Requirements according to BS EN 61000 and BS EN 60601-1:2006+A13:2024: Guidance and manufacturer's declaration – electromagnetic immunity			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s)	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT† (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
<i>†UT is the a.c. mains voltage prior to application of the test level.</i>			
Conducted RF IEC 61000- 4-6  Radiated RF IEC 61000- 4-3	3 Vrms 150 kHz to 80 MHz  3 V/m 80 MHz to 2.5 GHz	3 Vrms  3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the control unit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.  Recommended separation distance. d = 1.2VP d = 1.2VP 80 MHz to 800 MHz d = 2.3VP 800 MHz to 2.5 GHz

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

			<p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey*, should be less than the compliance level in each frequency range**.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
--	--	--	--

\* Field strengths from fixed transmitters cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Dynamic mattress system is used exceeds the applicable RF compliance level above, Dynamic mattress system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the system.

\*\* Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

<b>NOTE</b>	<p>At 80 MHz and 800 MHz, the higher frequency range applies. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>
-------------	---

Dynamic mattress system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer/user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Dynamic mattress system as recommended below, according to the maximum output power of the communications equipment.

<p><b>Recommended separation distances between portable and mobile RF communications equipment and the control unit</b></p>
---

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m) Electromagnetic environment – guidance		
	150 KHZ TO 80 MHZ D = 1.2√P	80 MHZ TO 800 MHZ D = 1.2√P	800 MHZ TO 2.5 GHZ D = 2.3√P
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23


For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE**

At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Pure Air 8 Auto Reactive (AR) Mattress System (PUR-8-AR)

### 17. Contact Information



WINNCARE PAC Ltd.  
Unit 2 Interchange Way,  
Oakenshaw  
Bradford,  
BD12 7AZ  
United Kingdom

Customer Service  
Tel +44 0800 193 6030  
Email: [askus@winncare.uk](mailto:askus@winncare.uk)  
Servicing Email: [askservice@winncare.uk](mailto:askservice@winncare.uk)  
[www.winncare.uk](http://www.winncare.uk)



WINNCARE France  
4, Le Pas du Chateau,  
85670 Saint-Paul-Mont-Penit  
France  
Tel. +33 (0)4 66 02 15 15

[www.winncare.fr](http://www.winncare.fr)



**NOTE**

Inform the competent authority if you believe or have reason to believe that the device presents a serious risk or that it has been tampered with.  
All serious incidents which are related to the device must be notified to the manufacturer and the competent authority of the member state in which the user and/or patient resides.