



MANGAR Leglifter

The inflatable Leglifter raises, lifts and transfers the users' legs into bed.

EASY TO FIT

- Fits most beds including profiling beds
- Straight forward to fit or remove, and requires no permanent installation
- Attaches to the bed underneath the mattress
- One airline connects the Leglifter to the Airflo Compressor and inflation is activated by pushing one button

SIMPLE TO MAINTAIN

- For a routine clean, sponge clean with a detergent solution and dry thoroughly
- · The covers can be machine washed
- · Spare covers can be purchased

 The materials are designed to exceptionally high standards with extraordinary durability, strength & performance in the toughest of environments

SAFE TO USE

- Raises, lifts and transfers legs into the bed
- Leg break on profiling beds can still be used once Leglifter is fitted
- · Aids independence

IMPROVES WELL-BEING

- Minimises moving and handling so reducing the risk of musculoskeletal injury for the carer
- Fits discreetly against the side of the bed
- · No permanent installation required

TECHNICAL SPECIFICATIONS

Max user weight	135 kg
Min height	Flat
Total length	124 cm
Width	57 cm
Total product weight	1.5 kg
Airflo Compressor	12
Lifting performance per charge	Up to 16 lifts

Features

- 1. Fitting Straps.
- 2. Air Distribution Hose.
- 3. Made from strong polyurethane coated nylon material.
- 4. Seams joined by high frequency welding.
- 5. Hook of Leglifter straps.







Product code: MPCA140500





For more information visit

www.mangarhealth.com Email: enquiries@winncare.uk Phone: 0800 2800 485



