


Product News


March 2021



WINNOCARE
N o r d i c

Smart Cart / Room-to-room spreader bar

SMART CART



UNIQUE XY-CART

WinnCare Nordics has developed a completely unique XY-cart, Smart Cart, which ensures that the healthy work environment already starts at the installation process.


Thanks to the Smart Cart, installation of a full room covering ceiling hoist system can now be done by only one installer.

This new solution is also time-saving, not only in installation context but also when it comes to order and expedition time.

Click on the brochure or the film for more details and information.



ROOM-TO-ROOM SPREADER BAR



A MULTI FUNCTIONAL SOLUTION

Optimized version of (and replacement for) WinnCare Nordics straight spreader bar.

The new spreader bar for Luna 275 kg has been designed with a unique 360° room-to-room system and can be changed without any use of tools.

Available in 360/480/580 mm
See item numbers below.



Click to see complete overview of WinnCare Nordic's spreader bars for Luna ceiling hoist

360 mm	480 mm	580 mm
015-02117	015-02118	015-02119

WINCARE
N o r d i c