

PITCHOUNE ROOM



DESIGNED WITH CHILDREN IN MIND...



PITCHOUNE Atmosphere
White lacquer surround, Kiwi or Violet inside finish
Room in standard colors

Chair not available for sale

FURNITURE FOR THE PITCHOUNE ROOM




BOARDS FOR COLLECTIVE CARE SETTINGS

- DIMENSIONS:
W953 x H553mm
- ▶ White lacquer surround
 - ▶ Kiwi or Violet melamine finish



BOARDS (FOR WOODEN SIDE RAIL)

- DIMENSIONS:
W932 x H636 mm
- ▶ White lacquer surround
 - ▶ Kiwi or Violet melamine-faced inside board or Orchid Violet



BEDSIDE CABINET

- DIMENSIONS:
W448 x D410 x H700 mm
- ▶ 1 drawer
 - ▶ Hinge door hung to the right
 - ▶ 1 shelf with bottle space




DRESSER

- DIMENSIONS:
W870 x D575 x H775 mm
- ▶ On feet
 - ▶ 3 drawers



TABLE WITH DRESSER

- DIMENSIONS:
W1340 x D555 x H739 mm
- ▶ 4 right-hand drawers



SINGLE-DOOR WARDROBE

- DIMENSIONS:
W670 x D594 x H1805 mm
- ▶ left-hand hinge
 - ▶ 1 shelf
 - ▶ 1 clothes rail

ACCESSORIES



"MY GIRL" STICKER BOARD

"MY BOY" STICKER BOARD

PITCHOUNE	FURNITURE AND BOARD COLORS		SURROUND, EDGING AND HAND RAIL COLORS	MATERIAL
	STANDARD COLORS	OTHER COLORS		
BOARDS	Kiwi inside or Orchid parma	Choice of melamine finishes from our color chart	Ral 9016	Surround in MDF and inside finish in melamine-faced MDF
FURNITURE	Standard production: Kiwi tops Ral 9016 handles Ral 9016 lacquer			

1-year warranty*. *Under normal conditions of use and excluding wearing parts - please see instructions of use

Photos not contractually binding

WINCARE Group - 4, Le Pas du Château - 85670 Saint-Paul-Mont-Penit - France
Tél : +33 (0)2 51 98 55 64 - Fax : +33 (0)2 51 98 59 07 - E-mail : info@medicatlantique.fr - www.winnicare.fr

In the interest of constantly improving its products, the WINCARE Group reserves the rights to modify the technical specifications of any product without prior notice
Class I medical device in accordance with European Directive 93/42/EEC

