



**Rise™ IIIS/130**  
Transfer, shower,  
commode and care chair  
with electrical elevation  
& seat tilt and soft  
backrest.

**Rise™ IIIS/130**  
Improved work  
environment and  
superb user comfort

# Rise™ IIIS/130

## Increased user comfort and better work safety

Powered, height adjustable, tilt-  
in-space shower commode chair



Now with  
37° seat tilt  
and 130 kg /  
20,7 st user  
weight

**BORRINGIA**

# Rise™ IIIS/130

## Increased user comfort and better work safety

### Powered, height adjustable, tilt-in-space shower commode chair

The Rise range of integrated care chairs set the standard, when the first model was introduced nearly 15 years ago. Based on years of user feedback, the Rise IIIS/130 introduces small but very important changes, making it the perfect tool in the care environment.

It is still incredibly compact and fits into even very small bathrooms. The electrical lift of 42 cm means that you can transfer the user directly from the bed without a person lifter

and guarantees the caregiver perfect working conditions throughout the care process.

In this way, straining of the caregiver's back, shoulders and knees is eliminated.

The lifting function allows the user to get up to the sink and mirror, where he can participate, thus creating a more dignified care routine.

With electrical seat tilt and multi-adjustable backrest, Rise IIIS /130 makes it easy for the user to find the most comfortable sitting position.



Rise™IIIS/130 with electric lifting function and tilt (backwards -30 / forward +7 degrees) Forward tilting seat in order to use the Rise to help the user standing up

### Key Features

- Electrical hoisting system (Soft Start/Stop)
- Electrical seat tilt (+7° / -30°) and multi-adjustable soft backrest offers increased user comfort and stability
- Removable battery w. external charging station
- Transfer directly from the bed without the use of a traditional person lift
- Fits over all standard toilets
- Suited for use in the shower
- Ideal working positions for the carer & increased efficiency
- Max User weight 130 kg / 286 lbs / 20,66 stone

# Optional Equipment – Rise™ III/130



Rise™ IIIS/130		Varenummer
	Rise IIIS/130	12C10S-130
Optional Equipment for Rise™ III/130		
1	Adjustable headrest, Small 12C16, Medium 12C17, Large 12C18	12C16/ 12C17/ 12C18
2	Headrest with flexible PU-padding	12C69
3	Swing-away Torso Support with PU padding	12C65
4	PU-toilet seat with narrow hole and lid, 12C22, Soft cover for the toilet seat 90D91	12C22/90D91
5	Calf strap for the footrests	90H27
6	125 mm wheels with brakes (1 pcs.) Can be used instead of the standard 100 mm wheels	90H86
7	Elevating legrests1) Set, 90C60 - 2) Left, 12C66 - 3) Right, 12C67	90C60/12C66/12C67
8	Hand control to RiseIIIS & Rise IIIS/130	90J09
9	Wide armrests (set). Increases the distance between the armrests by 10 cm	12C84
10	Hemi-armrests 1) Left 12C59 - 2) Right, 12C60	12C59/12C60
11	Locks on the armrests (set)	12C68
12	Amputation Leg rest, 1) Left, 12C63- 2) Right, 12C64	12C63/12C64
13	Positioning backrest with built-in pillow	12C52
14	Bracket for wall mounting of charging station	12C48
15	Seat Tilt Limitation. Blocks tilt forward	12C55
16	Crossbar between the armrests. Padded. Secured with elastic bands	12C61
17	Universal Anti- Pressure Pillow for toilet seat. Pillow & covers can be disassembled, washed at 90 degrees. And tumble-dried	90A19S
18	Bedpan Holder for low hanging bedpan. Extends the distance to the seat by 50 mm	12C68
19	Bedpan Holder with stops that hold the bedpan in place, when the seat is tilted.	12C53
20	Direction lock (factory mounted)	12C75
21	Torso-, Hip & Leg Belt	90E08
22	Extra Battery 90J08, Extra Charging Station 90G77	90J08/90G77
23	Bedpan with lid, stainless steel, 12A09, Bedpan without lid, 12C56, Lid for bedpan 12C57	12A09/12C56/12C57

Rise™ IIIS/130 is CE-Marked according to ISO/EN 10535. Rise™ IIIS/130 is patented. Rise IIIS/130 is a trademark ™, belonging to the manufacturer.



## Technical specifications

<b>Max. capacity (SWL)</b>	130 kg / 20,7 st	<b>Lifting machinery</b>	Electrical Actuator
<b>Height</b>	1065 - 1480 mm	<b>Battery</b>	24 volt
<b>Width</b>	570 mm	<b>Battery capacity</b>	25-30 løftecykler
<b>Length</b>	880 - 1050 mm	<b>Charging</b>	External charging station
<b>Seat height</b>	535 - 950 mm	<b>Charging time</b>	4-5 hours
<b>Seat width</b>	450 mm	<b>Turning circle</b>	980 mm
<b>Seat depth</b>	430 mm	<b>Lowest clearing height</b>	70 mm
<b>Weight</b>	32 kg	<b>Chassis height</b>	150 mm



MDD93/42/EEC, Class I

Consult the User manual for further specifications. All figures and measurements are for guidance only and may differ in actual situations. Winncare reserves the right to changes without notice. Product photos may include accessories. Consult your dealer.

## "From-2-to-1-Care Givers"

With the increasing number of high dependency elderly people, we will soon be missing warm hands, if we keep working in the traditional manner. Luckily, a number of independent studies in Denmark have shown that it is possible to go "From-2-to-1-Care Givers" in most cases.

By combining the Luna Ceiling Lifter with the Rise chair, one caregiver can safely handle all

routines in connection with lifting/transfer from the bed, toileting, showering, general hygiene and getting dressed. The experience from a number of Danish municipalities show

- Better, more dignified care
- Less stress, better work environment and reduced sick leave
- Increased efficiency and large cost savings

**Please contact us, if you want to learn more about "From-2-to-1-Care Givers".**



by Winncare Nordic

Winncare Nordic ApS · Taarnborgvej 12C · 4220 Korsoer · Denmark  
Tel: +45 70 27 37 20 · e-mail: [ergolet@ergolet.dk](mailto:ergolet@ergolet.dk)

[www.ergolet.com](http://www.ergolet.com)