



## MANGAR Evac

The Mangar Evac is a flexible, compact stretcher system for use by the emergency services, to evacuate immobile people in the event of a medical emergency or other emergency situation where the use of conventional patient transfer equipment would not work.

### MINIMISING IMPACT

The Mangar Evac holds the patient securely but with the flexibility to negotiate obstacles. It is designed to facilitate the manual movement of a patient by maximising the use of sliding techniques where possible and minimising the impact on the person performing the movement.

### SEAMLESS TRANSFER

Where required, the Evac can also be secured around the patient whilst still on the Rhino and/or Swift. This reduces the number of transfers required to move a patient from a bed into the Evac.

### ULTIMATE VERSATILITY

The Mangar Evac can navigate obstacles, narrow hallways, slopes, and stairways, and has been used by emergency crews in many different scenarios to achieve a safe transfer, including being hoisted by a helicopter.



Mangar Rhino Lifting Cushion  
and Swift Transfer Slide



SCAN ME



## SPECIFICATION

- **Packed size:** 70cm x 35cm x 35cm
- **Weight:** 14kg
- **Deployed size:** 230cm x 260cm x 2cm
- **Safety tested to:** 1000kg
- **Working Load Limit (WLL):** 500kg
- **Material Construction:** P300 Plastic base (black 2mm board is replaceable), heavy duty PVC upper, 50mm colour coded safety webbings fitted with metal seatbelt type adjustable buckles.

## KEY FEATURES

- Removeable base can be replaced if damaged
- Colour coded webbed straps
- Anchor loops at head and foot end can be used for hoisting or winching patients
- Compatible with the Mangar Rhino and Swift for a seamless transfer
- Flame-retardant PVC material can be wiped or washed clean after each use
- Manufactured in the UK



**For more information visit**

[www.winnccare.uk](http://www.winnccare.uk)

Email: [enquiries@winnccare.uk](mailto:enquiries@winnccare.uk)

Phone: 0800 2800 485