



LUNA[®] X5


Accessories



Spreader bars



2-point



Small

Size: 370 mm
Item number: 017-00335




Medium

Size: 470 mm
Item number: 017-00336



Large

Size: 570 mm
Item number: 017-00337



Holder for 2-point spreader bars (size S-M-L)

L = 235 mm H = 47 mm
Item number: 016-02230

Room-to-room



For room-to-room transfers using two LUNA® X5

Item numbers:
360 mm: 015-02117
480 mm: 015-02118
580 mm: 015-02119

Use with bayonet adapter - Item number: 017-00331

3-point



Enhanced patient comfort

For larger clients, WinnCare's V-shaped 3-point spreader bar increases the distance between the shoulder and leg sections, opening the sling for a more comfortable lift.

Item number: 017-00338

4-point



Support for bariatric transfers

For bariatric transfers, we recommend our 4-point spreader bar. Combined with the correct sling model and size, it provides proper support and comfort for larger users.


Item number: 017-00342

3-Point Powered Positioning spreader bar



DIMENSIONS

POSITION:	SITTING:	LYING:
A) Track to sling attachment point:	935 mm	725 mm
B) Height:	695 mm	475 mm
C) Width:	120 mm	520 mm



Enhanced patient comfort

- SWL 275 kg
- 360 degree rotation
- Easy attachment to LUNA® X5
- Simple and intuitive positioning via the LUNA® X5 hand control

Item number: 017-00340

Use with Winncare's Click Clip Sling



Trolley



Lightweight trolley for moving LUNA® X5 from one room to another

- On board charging
- Holder for extension arm

Item number: 017-00354

Weighing scale



- **Medically approved weighing scale**
The weighing scale is designed to weigh disabled patients, when they are placed in a sling.
 - Battery powered (AAA)
 - On/Off/Zero/Tare/Hold button
 - Weighing capability up to 300 kg
- Item number: 017-01422

Bayonet adapter



- **Adapter for bayonet hook**
Use with room-to-room spreader bar or for attaching the IPX4 weighing scale to the LUNA® X5.
- Item number: 017-00331




Extension strap



- **Option for systems intalled higher than three meters**
Item numbers:
500 mm: 017-00154
1000 mm: 017-00153

Extension arm



- **For installing LUNA® X5 in portable mode**
 - Telescopic extension arm to help mounting the LUNA® X5
 - Available in two sizes:
650 mm - 1100 mm: Item number: 017-00350
1000 mm - 1800 mm: Item number: 017-00351

Hand control / Side panel



• User friendly hand control (IPX7)

Available in standard and two channel.
Use the two channel hand control for operating the 3-point Powered Positioning spreader bar.

Item numbers:

Standard + 1080 mm cable: 017-00320

Standard + 1320 mm cable: 017-00321

Two channel + 1080 mm cable: 017-00322

Two channel + 1320 mm cable: 017-00323



• Two channel side panel kit

Use the two channel side panel for operating the 3-point Powered Positioning spreader bar.


Item number: 017-00330

Charging



• Charging station

- 1) Charging station: 017-00301
- 2) Charger with plug and 1080 mm cable: 017-00302
- 3) Wall bracket set for charging station: 017-00303
- 1+2+3: Combined charging set: 017-00300
- 4) Extension cable for charger, 1850 mm: 017-00304



• In-Track Charging (ITC) of the LUNA® X5

Consists of the following item numbers:

ITC plastic track: 016-01207

ITC tin-bronze strip: 016-01206

ITC conductor / end cover: 016-01199

ITC conductor / LUNA® X5 cart: 016-01201

ITC conductor / XY cart: 016-01202

ITC plastic track ramp: 016-01208

Charger for NiMH batt./1,8m DC cable: 017-00302xx*

*EU / UK / AUS / US



winncare.com