



DIGNIFIED CARE



IMPROVED WORKING ENVIRONMENT



EFFICIENT AND TIME-SAVING

TRITON DUAL 150

Powered, height adjustable, tilt-in-space shower commode chair

BETTER WORKING ENVIRONMENT AND MORE COMFORT

- Correct working posture eliminates strain and stress on the caregiver's back, shoulders and knees.
- With the powered seat tilt and adjustable backrest, Triton Dual 150 makes it easy for the user to find the most comfortable sitting position.
- The height adjustment allows an optimal positioning of the chair, ensuring a correct working posture for the caregiver.
- For the user, the height adjustment means increased independence; it is possible to adjust the chair, so the user can reach the sink and look in the mirror - thereby the user can help take care of his own personal hygiene.

FEATURES

- If the seat is raised while tilted forwards, the Dual Security function will simultaneously tilt the seat backwards (5 degrees). This way, the user will not have the feeling of falling out of the chair while it is moving.
- Transport function: Press down the transport button on the hand control, and the chair is positioned for safe transport of the user.
- Triton Dual 150 also comes with two memory functions which makes it quick and easy for the carer to adjust the chair in a desired position, for example for hair wash or toilet visits.

TECHNICAL SPECIFICATION

Max. user weight: 150 kg
(330 lbs/23,62 stone)

Distance btw. arm rest: 470 mm

Seat depth: 430 mm

Seat width: 460 mm

Seat height: 552-930 mm

Total width: 620 mm

Length: 840-1000 mm

Total height: 1119-1497 mm (without head rest)

Weight: 38 kg

Tilt angle: +7° til -30° ±2°

IP class: IPx6W

Charging time: 4 hours

Battery capacity: 80-100 lifting cycles

PRODUCT HIGHLIGHTS

- Powered height adjustment (Soft Start/Stop).
- Powered tilt (+7° / -30°) and adjustable backrest for more comfort and user stability.
- IPx6W rated: Can be cleaned in washing machines up to 50 times due to LINAK IPX6 Washable™
- Developed from the Easy-Clean principle, which means there are no small corners or edges. The chair has closed screws/joinings and fully welded steel, which prevents water from entering the pipes. The lifting actuator is covered, making it easy to clean in case of splashes or impurities.
- Treated with double coating which makes it rust resistant.
- Dual Security Function (it can not be raised and tilted forward simultaneously)
- Transport function via the hand control for safe transportation of the user.
- Two programmable memory buttons on the hand control for automatic positioning and adjustment.
- Multi adjustable headrest with integrated rod which does not scratch or damage the wall behind the toilet
- Fits over all standard toilets and suitable for use in the shower
- The battery must be charged regularly with the attached cable or in an external charging station.
- Max user weight 150 kg



Please consult the technical manual for further specifications. All figures and measurements are for guidance only and may differ in actual situations. WinnCare Nordic reserves the right to changes without notice. Product photos may include accessories. Please consult your dealer for further information.

FOR MORE INFORMATION VISIT

winnCare.com
E-mail: info@winnCare.dk
Telefon: +45 70273720

