

ALDRYS'EVOLUTION II[®]



Bed with CARMEN II boards. Light Dragon Tree finish, RAL1013 handrail.



AN ELDERLY PERSON'S NEEDS CAN EVOLVE... SO CAN OUR BED

STEP N°1: Choose the bed base according to the variable height



80 cm
24 cm
if centralized braking
Ø 100 mm



STEP N°2: Choose the bed configuration

▶ BACK REST: several options available



Auto-regression back rest.
It prevents the patient from sliding forward.



Emergency release for the back rest
can be accessed from the head and
both sides of the bed.

▶ LEG REST: several options available



Electric leg rest with thigh section.



▶ BRAKING



Castors Ø 75 mm
Independent locking
on each castor.

Locked position



Unlocked position



The single pedal is located at the foot
of the bed in the middle.

HEIGHT RANGE

- ▶ 22 to 78 cm
- ▶ Safe for patients, comfortable for caregivers.

OPTIMIZED UNDER BED CLEARANCE

ROBUST

- ▶ LINAK motorization and materials: 5-year warranty.

SAFETY

- ▶ Colored, traditional or padded (Elegantti) bedends to give patients a visual cue and extra safety in the event of a bump.

DESCRIPTION	DIMENSIONS	ADJUSTABLE HEIGHT	BACK REST	LEG REST	LIFTING CAPACITY	STANDARDS	MOTORIZATION
ALDRYS'EVOLUTION II 2 electrical functions 3 electrical functions without knee break	90 x 200 cm divisible bed base or non-divisible	22 CM à 78 cm freinage roue individuelle 24 cm à 80 cm freinage centralisé	70°	12° manual cranks	SWL ⁽¹⁾ 170 kg mattress 20 kg patient 135 kg accessories 15 kg	IEC 60601-2-52 compatible with environments 3, 4 and 5	LINAK IP54
3 electrical functions with thigh sections				Knee break 26° (thighs) and 16° feet adjustment by manual cranks			

⁽¹⁾ Safe working load

- 2 accessory holders (lengthwise)
- 2 IV stand holders (at the head of the bed)
- 2 lifting pole (at the head of the bed)

Handset, electric and mechanical parts: 5-year warranty*

Boards and accessories: 1-year warranty*

* Under normal conditions of use and excluding wearing parts - please see instructions of use

Photos not contractually binding

FP000003/01.2018

